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
Analytes	LEGENDplex™ Sensitivity (pg/mL)*	ELISA Sensitivity (pg/mL)**
IFN- γ	1.1	4
TNF- α	0.7	3.5
IL-2	1.4	4

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LEGENDplex™ Mouse Anti-Virus Response Panel (13-plex) with Filter Plate	Y	Multiplex - Panel	RUO	740621	100 tests
LEGENDplex™ Mouse Anti-Virus Response Panel (13-plex) with V-bottom Plate	Y	Multiplex - Panel	RUO	740622	100 tests
LEGENDplex™ Mouse Cytokine Panel 2 (13-plex) with Filter Plate	Y	Multiplex - Panel	RUO	740134	100 tests
LEGENDplex™ Mouse HSC Erythroid Panel (7-plex) with Filter Plate	Y	Multiplex - Panel	RUO	740680	100 tests
LEGENDplex™ Mouse HSC Erythroid Panel (7-plex) with V-bottom Plate	Y	Multiplex - Panel	RUO	740681	100 tests
LEGENDplex™ Mouse HSC Lymphoid Panel (7-plex) with Filter Plate	Y	Multiplex - Panel	RUO	740678	100 tests
LEGENDplex™ Mouse HSC Lymphoid Panel (7-plex) with V-bottom Plate	Y	Multiplex - Panel	RUO	740679	100 tests
LEGENDplex™ Mouse HSC Myeloid Panel (7-plex) with Filter Plate	Y	Multiplex - Panel	RUO	740682	100 tests
LEGENDplex™ Mouse HSC Myeloid Panel (7-plex) with V-bottom Plate	Y	Multiplex - Panel	RUO	740683	100 tests
LEGENDplex™ Mouse HSC Panel (13-plex) with Filter Plate	Y	Multiplex - Panel	RUO	740676	100 tests
LEGENDplex™ Mouse HSC Panel (13-plex) with V-bottom Plate	Y	Multiplex - Panel	RUO	740677	100 tests
LEGENDplex™ Mouse Inflammation Panel (13-plex) with Filter Plate	Y	Multiplex - Panel	RUO	740150	100 tests

Humain

LEGENDplex™ HU Essential Immune Response Panel (13-plex) w/FP

Pricing & Availability

[Kit Contents](#)

[Product Details](#)

[Antigen Details](#)

[Documentation](#)

[Reviews](#)

Regulatory Status RUO

Other Names IL-4, IL-2, CXCL10 (IP-10), IL-1β, TNF-α, CCL2 (MCP-1), IL-17A, IL-6, IL-10, IFN-γ, IL-12p70, CXCL8 (IL-8), TGF-β1

Ave. Rating [Submit a Review](#)

Product Citations [1 publications](#)

Souris

LEGENDplex™ Mouse Inflammation Panel (13-plex) with Filter Plate

Pricing & Availability

[Kit Contents](#)

[Product Details](#)

[Antigen Details](#)

Regulatory Status RUO

Other Names IL-1a, IL-1b, IL-6, IL-10, IL-12 p70, IL-23, IL-27, IFN-b, IL-17A, MCP-1, TNF-a, IFN-g, GM-CSF, IL-1 alpha, IL-1 beta, IFN beta, TNF alpha, IFN gamma, il1 alpha, il-1 alpha, il1alpha, il1 a, il-1 a, il1a

Rat

LEGENDplex™ Rat Inflammation Panel (13-plex) with Filter Plate

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
Overview	Panel Wizard	Pre-defined Panels	Mix and Match	Custom Panels
Software	Instrument Setup	Hands-on Training	Support	Ancillary Products

Click on the buttons below to learn more about our LEGENDplex software solutions.

[Data Analysis Software Suite](#)[Legacy Programs](#)

The LEGENDplex™ Data Analysis Software Suite is designed to analyze LEGENDplex™ flow cytometry data files.

The software suite consists of both an online cloud-based program as well as downloadable desktop applications for both PC and Mac computers.

Please click the User Manual button below for full details on software suite usage and installation.

Shared Core Software Features:


- A simple, intuitive user interface and workflow
- Pre-loaded LEGENDplex™ assay bead IDs
- A sophisticated auto-gating algorithm
- Simple import of experiment workspaces between cloud-based program and desktop applications

Cloud-Specific Features:

- Access from any web-connected computer
- 2 GB of storage space per user account
- One click access to previously uploaded data
- Easily share analysis results with other users

[Access the software suite in a new browser window](#)

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The screenshot shows the LEGENDplex website. A black arrow points to the 'Software' tab in the top navigation bar. Below the navigation bar, there is a section titled 'Click on the buttons below to learn more about our LEGENDplex software solutions.' This section contains two buttons: 'Data Analysis Software Suite' and 'Legacy Programs'. Below these buttons is a video player showing a screenshot of the LEGENDplex Cloud-Based Software Suite interface. Below the video player is a large play button icon. Below the play button is a screenshot of the LEGENDplex software suite running on a desktop. Below the screenshot is a paragraph of text: 'The LEGENDplex™ Data Analysis Software Suite is designed to analyze LEGENDplex™ flow cytometry data files. The software suite consists of both an online cloud-based program as well as downloadable desktop applications for both PC and Mac computers. Please click the User Manual button below for full details on software suite usage and installation.' Below this text are two columns of features: 'Shared Core Software Features' and 'Cloud-Specific Features'. At the bottom of the page is a button that says 'Access the software suite in a new browser window', with a black arrow pointing to it from the right.

Overview	Panel Wizard	Pre-defined Panels	Mix and Match	Custom Panels
Software	Instrument Setup	Hands-on Training	Support	Ancillary Products

Click on the buttons below to learn more about our LEGENDplex software solutions.

Data Analysis Software Suite Legacy Programs

LEGENDplex™ Cloud-Based Software...

The LEGENDplex™ Data Analysis Software Suite is designed to analyze LEGENDplex™ flow cytometry data files.

The software suite consists of both an online cloud-based program as well as downloadable desktop applications for both PC and Mac computers.

Please click the User Manual button below for full details on software suite usage and installation.


Shared Core Software Features:

- A simple, intuitive user interface and workflow
- Pre-loaded LEGENDplex™ assay bead IDs
- A sophisticated auto-gating algorithm
- Simple import of experiment workspaces between cloud-based program and desktop applications

Cloud-Specific Features:

- Access from any web-connected computer
- 2 GB of storage space per user account
- One click access to previously uploaded data
- Easily share analysis results with other users

Access the software suite in a new browser window

The screenshot shows the login form on the LEGENDplex website. It has a 'Log In' header. Below the header are two input fields: 'Email' and 'Password'. Below the 'Password' field is a link that says 'Forgot Password'. At the bottom of the form is a large purple button that says 'Login'.

Log In

Email


Password

[Forgot Password](#)

Login

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Inflammation Panel (13-plex) (Mouse)

Inflammation Panel (13-plex) (Mouse)

Inflammation Panel (13-plex) (Rat)

Inflammation Panel 1 (13-plex) (Human)

Inflammation Panel 2 - S/P (11-plex) (Human)

Inflammation Panel 2 - S/P (2-plex) (Human)

Inflammation Panel 2 - TC (13-plex) (Human)

Kidney Function Panel 1 (9-plex) (Human)

Kidney Function Panel 2 (7-plex) (Human)

Kidney Function Panel 3 (5-plex) (Human)

M1 Macrophage Panel (8-plex) (Human)

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TOP STANDARD CONCENTRATIONS AFTER RECONSTITUTION:

Analyte	Bead ID	Concentrations (ng/mL)	Analyte	Bead ID	Concentrations (ng/mL)
IL-23	A4	96	IL-1 β	B2	12
IL-1 α	A5	15	IL-10	B3	15
IFN- γ	A6	28	IL-6	B4	16
TNF- α	A7	15	IL-27	B5	72
MCP-1	A8	14	IL-17A	B6	16
IL-12p70	A10	16	IFN- β	B7	300
			GM-CSF	B9	15

This lot of product has met all specifications for release.

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The screenshot shows the 'Samples' tab of a software interface. At the top, there are tabs for 'Beads', 'Standards', 'Samples' (selected), 'Review Gating', and 'Results'. Below the 'Samples' tab, there is a 'Replicate Controls' panel. Inside this panel, the 'Replicates' field is set to 2, and the 'Replicate Pattern' is set to 'AABBC'. An 'Apply' button is located below these fields. Below the 'Replicate Controls' panel, there is a button labeled 'Select Samples FCS file(s)' with a folder icon. An arrow points from the 'Ajout des échantillons' header to this button. Below the 'Select Samples FCS file(s)' button, there is an 'FCS File List' section. It includes a table with the following columns: 'Files', 'Sample Name', 'Dilution', and 'SBTIM (Acquisition Time)'. The table is currently empty. Below the table, there is an 'Upload' button.

Files	Sample Name	Dilution	SBTIM (Acquisition Time)
-------	-------------	----------	--------------------------

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Replicate Pattern:

FCS File List:

Files	Sample Name	Dilution	SBTIM (Acquisition Time)

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


Uploaded Sample and Standard Files:


Files	Sample Name	Replicate	Dilution	
S101-2.fcs	S101-2	1	2	✓
S102-1.fcs	S102-1	1	2	✓
S102-2.fcs	S102-2	1	2	✓
S103-1.fcs	S103-1	1	2	✓
S103-2.fcs	S103-2	1	2	✓
S104-1.fcs	S104-1	1	2	✓
S104-2.fcs	S104-2	1	2	✓
S105-1.fcs	S105-1	1	2	✓
S105-2.fcs	S105-2	1	2	✓
S106-1.fcs	S106-1	1	2	✓
S106-2.fcs	S106-2	1	2	✓
S107-1.fcs	S107-1	1	2	✓
S107-2.fcs	S107-2	1	2	✓

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
Color By: ☐ Concentration ☒ Count ☐ MFI ☐ CV

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